# VISALYS® FILL

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## INSTRUCTIONS FOR USE USA | EN

In accordance with DIN EN ISO 4049  
Carefully read instructions prior to use

### Product description

Visalys® Fill is a light-cured filling material combining proven composite and innovative nano technology. Visalys® Fill is radiopaque, contains 87 % w/w inorganic fillers (= 71.4 Vol. %) in a methacrylate matrix (BIS-GMA, TEDMA, UDMA) and cures under blue light (halogen or LED).

Visalys® Fill is used with a dentin/enamel bond. Visalys® Fill is available in syringes and in caps for direct application.

### Shades

A1, A2, A3, A3.5, A4, OA2, BL

### Indications

- Filling class I to V cavities
- Reconstruction of traumatically affected anteriors
- Facetting of discolored anteriors
- Correcting shape and shade for better esthetic appearance
- Locking, splinting of loose anteriors
- Repairing veneers
- Filling deciduous teeth
- Core build-up under crowns
- Composite inlays

### Contraindications

Visalys® Fill should not be used if the recommended application technique cannot be used or if the patient has a known allergy to any of the ingredients (BIS-GMA, BHT, camphor quinone, UDMA, TEDMA).

### Preparation / Shade selection

Clean the teeth to be treated. Mark occlusal contact points (posterior area); a minimal separation facilitates the design of the approximal contact and placing of the matrix. Clean the teeth prior to shade selection. The shade is selected in comparison with the tooth while it is still moist.

### Drying

Ensure that the working field is dry. Using a rubber dam is recommended.

### Cavity preparation

Generally, cavity preparation should be carried out according to the rules of the adhesive filling therapy, i.e. minimally invasive to conserve healthy tooth substance. Bevel the enamel margins on anteriors and round off the preparation margins on posteriors. Non-carious cervical lesions do not have to be prepared; thorough cleaning is sufficient here. Afterwards, clean and dry the cavity.

### Lining

A calcium hydroxide lining should be applied in proximity of the pulp. Place a layer of stable cement over the lining. Due to their fluoride release, glass polyalkenoate (glass ionomer) materials are recommended.

### Bond material

Visalys® Fill is used in adhesive technique with a dentin/enamel bond. All light-curing bonding materials may be used. Follow the respective instructions for use with regards to preparation (etch technique) and application.

### Application of Visalys® Fill

Let the material reach room temperature before application.

Place a matrix on cavities in the approximal area. Using translucent matrices is advantageous (smooth surfaces). Apply the chosen shade of Visalys® Fill (see shade selection) and adapt it with a suitable instrument. Apply and polymerize fillings of more than 2 mm in layers.

### Light-curing

Conventional polymerisation devices are suited for light-curing this material. The polymerisation time per layer (Polymerisation time) as a function of the light intensity (Light intensity) can be found in Table 1.

The curing time is a minimum of 20 sec per layer, for opaque shade OA2 a minimum of 40 sec. Hold the light emission tip of the device as close as possible to the surface of the filling. If the distance is more than 5 mm, the curing depth may be compromised. Incomplete curing may lead to discoloration and pulpitis-like complaints.

Table 1

Light intensity (LED-/halogen lamp)	Polymerisation time (per layer)	
≥ 500 mW/cm <sup>2</sup>	20 s A1, A2, A3, A3.5, A4, BL	40 s OA2
≥ 800 mW/cm <sup>2</sup>	10 s A1, A2, A3, A3.5, A4, BL	20 s OA2

### Finishing

The filling can be finished and polished immediately after removing the shaping aids (e.g. fine or extra fine diamond burst, polishing disc), with cooling provided. The margin of the filling or the entire tooth should be fluoridated as a final step.

Visalys® Fill can be used for direct and indirect inlays according to customary methods. The physical stability can be improved by the usual external (additional) curing.

### Information / Precautionary measures

- No known side effects. Hypersensitive persons may develop sensitivities.
- Phenolic substances, especially eugenol and thymol-containing preparations, interfere with curing filling composites. The use of zinc oxide-eugenol cements or other eugenol-containing materials in combination with filling composites should be avoided.
- Fillings exposed to occlusal forces should be checked at least once a year for early detection of changes.

### Storage

Store at temperatures between 39°F – 82°F (4°C – 28°C). If refrigerated, the material must be allowed to reach room temperature before use. To avoid exposure to light and possible polymerization, syringes should be closed immediately after dispensing. Do not use Visalys® Fill after the date of expiration.

This material has been developed solely for use in dentistry. Processing should be done strictly according to the instructions for use. Kettenbach recognizes its responsibility to replace products if proven to be defective. Kettenbach does not accept liability for any damage or loss, directly or indirectly, stemming from the use of or inability to use the products described. Before using, it is the responsibility of the user to determine the suitability of the product for its intended use. The user assumes all risk and liability in connection therewith.

Descriptions and data constitute no warranty of attributes and are not binding.

### CAUTION: U.S. Federal Laws restrict this device to sale by or on the order of a dentist.

No person is authorized to provide any information which deviates from the information provided in the instructions for use.

Keep this material out of reach of children.

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Fabriqué pour : Kettenbach GmbH & Co. KG