

PREVENTION

DIGITAL
RESTORATION

RESTORATION

TMD

IMPRESSIONS

FUTAR® – BITE REGISTRATIONS MADE EASY

A-SILICONE FOR BITE REGISTRATIONS

UNIQUE: 100% DEVELOPED AND MANUFACTURED IN GERMANY

We oversee every step – from initial concept to final product!
True German manufacturing, true responsibility, and quality you can truly trust.



MADE IN
GERMANY



www.kettenbach-dental.us

Strong in prevention – restoration – impressions.


KETTENBACHDENTAL
Simply intelligent

PROVEN MATERIALS. CONVENIENT AND PRECISE.

FUTAR® D SLOW

- Extended working time (90 seconds) for myocentric bite registration.
- Good machinability.

FUTAR® D / FUTAR® D FAST

- The most widely used product in Germany!
- Excellent millability and precise results guaranteed – thanks to the high final hardness!

FUTAR® / FUTAR® FAST

- On the market since 1994, making it the classic of the product family.
- Hard after setting, flexible in processing: very good for trimming.

FUTAR® EASY / FUTAR® EASY FAST

- Hard material with low dispensing force – easy to apply.
- Good for trimming.

FUTAR® CLEAR / FUTAR® CLEAR FAST

- Direct bite control thanks to transparent material.
- Convenient final hardness of A 62.

Technical data

TOTAL SETTING TIME*

1:00 MIN. –
SEGMENT

2:00 MIN. –
DENTAL ARCH

3:00 MIN. –
DENTAL ARCH

4:30 MIN. –
DENTAL ARCH



1:1 cartridge

FUTAR® FAST (Shore A 90)
FUTAR® D FAST (Shore A 43)
FUTAR® EASY FAST (Shore A 90)
FUTAR® CLEAR FAST (Shore A 62)

FUTAR® (Shore A 90)
FUTAR® EASY (Shore A 90)
FUTAR® D (Shore A 43)

FUTAR® CLEAR (Shore A 62)

FUTAR® D SLOW (Shore D 43)

*Total setting time (removal from mouth) from start of mixing.

0211257_1926

Kettenbach LP

16052 Beach Blvd. Ste. 221 · Huntington Beach, CA 92647, USA · Phone: 877-532-2123 · Fax: 888-892-9820
Website: www.kettenbach-dental.us · Web shop: www.kettenbachusa.com · Custservice@kettenbachusa.com
Strong in prevention – restoration – impressions.



KETTENBACHDENTAL
Simply intelligent