

PREVENTION

DIGITAL
RESTORATION

RESTORATION

TMD

IMPRESSIONS

FUTAR® CLEAR & FUTAR® CLEAR FAST – BITE REGISTRATIONS MADE EASY

A-SILICONE FOR BITE REGISTRATIONS

**FUTAR CLEAR IS A TRANSPARENT, ELASTOMERIC
BITE REGISTRATION MATERIAL BASED ON A-SILICONE FOR
THE PRECISE REGISTRATION OF OCCLUSION.**

Its transparency allows direct visual control of the bite position during application.

Available in Regular and Fast versions – the Fast variant with reduced working time is particularly suitable for smaller applications and quick registrations.



**MADE IN
GERMANY**



www.kettenbach-dental.us

Strong in prevention – restoration – impressions.


KETTENBACHDENTAL
Simply intelligent

PROVEN MATERIALS. CONVENIENT AND PRECISE.

FUTAR® CLEAR

- Direct bite control thanks to transparent material
- No risk of fracture caused by comfortable final hardness (Shore-A 62)
- High thixotropy and flowability, no bite resistance
- Does not flow into the interdental spaces; easy removal from the mouth

FUTAR® CLEAR FAST

- Extra quick working time of 15 seconds, extra short intraoral setting time of 45 seconds
- Minimal time required, reduced risk of distortion during intraoral setting; more comfortable for the patient



YOUR FAVORITE
FUTAR® CLEAR FAST

NOW ALSO AS
FUTAR® CLEAR



Technical data



1:1 cartridge

TOTAL SETTING TIME*

1:00 MIN. – SEGMENT

3:00 MIN. – DENTAL ARCH

FUTAR® CLEAR FAST (Shore A 62)

FUTAR® CLEAR (Shore A 62)

*Total setting time (removal from mouth) from start of mixing.

0211257_1926

Kettenbach LP

16052 Beach Blvd. Ste. 221 · Huntington Beach, CA 92647, USA · Phone: 877-532-2123 · Fax: 888-892-9820
Website: www.kettenbach-dental.us · Web shop: www.kettenbachusa.com · Custservice@kettenbachusa.com
Strong in prevention – restoration – impressions.

KETTENBACHDENTAL
Simply intelligent