



 **SUPER WITH FUTAR®**

**BITE REGISTRATION**

  
**KETTENBACHDENTAL**  
Simply intelligent



# FUTAR® – SIMPLY IMPRESSIVE!

## Bite registration

The **Futar®** product family includes eight different bite registration materials based on A-silicone for making accurate occlusal records.

### **SUPER** IN QUALITY AND PROGRESS

- "Made in Germany" with its own research and development as well as production in Germany

### **SUPER** THROUGH VERSATILITY

- Choose between high and highest final hardness in addition to different processing times!

### **SUPER** FOR OVER 20 YEARS

- Trust Germany's market leader<sup>1)</sup>!

**Do not make compromises with your bite registration!**



<sup>1)</sup> Acc. to GfK, a German consumer research association

# FUTAR® – THE ADVANTAGES AT A GLANCE

NEW

## Futar® / Futar® Fast:

- High final hardness (Shore-A 90)
- Flexible in processing: easy to cut
- Also available with fast-setting characteristics: for more patient comfort!



## Futar® Easy / Futar® Easy Fast:

- High final hardness (Shore-A 90)
- Low bite resistance
- Easy to cut
- Also available with fast-setting characteristics
- Low dispensing force – easy to apply



The Futar® family can  
do even more:



## Futar® D / Futar® D Fast:

- Extra high final hardness (Shore-D 43)
- Offers very good millability
- Also available with fast-setting characteristics: for more patient comfort!



## Futar® D Slow:

- Extra high final hardness (Shore-D 43)
- Extra long processing time for time-consuming applications
- Versatile use also in the laboratory



- Accurate results thanks to excellent flow properties
- Low elasticity reduces springiness in the model assignment
- Surplus can be broken off easily









<sup>2)</sup> Total setting time (removal time from mouth) from the start of mixing

**SIMPLE,  
CLEAN AND  
FAST TO USE!**

# FUTAR® – PERFORMANCE RUNS IN THE FAMILY










## Dental arch

 Bite registration material	 Mixing				 Hardness (approx.) Shore
		Working time at 23 °C/74 °F ≤	Intraoral setting time at 35 °C/95 °F ≥	Total setting time <sup>2)</sup> ≥	
<b>Futar®</b>		30 sec.	90 sec.	<b>2:00 min.</b>	A 90
<b>Futar® Easy</b>		30 sec.	90 sec.	<b>2:00 min.</b>	A 90
<b>Futar® D</b>		30 sec.	90 sec.	<b>2:00 min.</b>	D 43
<b>Futar® D Slow</b>		90 sec.	3:00 min.	<b>4:30 min.</b>	D 43

<sup>2)</sup>Total setting time (removal time from mouth) from the start of mixing

## Segment

 Bite registration material	 Mixing				 Hardness (approx.) Shore
		Working time at 23 °C/74 °F ≤	Intraoral setting time at 35 °C/95 °F ≥	Total setting time <sup>2)</sup> ≥	
<b>Futar® Fast</b>		15 sec.	45 sec.	<b>1:00 min.</b>	A 90
<b>Futar® Easy Fast</b>		15 sec.	45 sec.	<b>1:00 min.</b>	A 90
<b>Futar® D Fast</b>		15 sec.	45 sec.	<b>1:00 min.</b>	D 43

<sup>2)</sup>Total setting time (removal time from mouth) from the start of mixing

CREATE  
YOUR BITE  
REGISTRATION  
IN JUST ONE  
MINUTE!



# FUTAR® – THE RANGE



Products	Normal pack 2 x 50 mL cartridges, 6 blue mixing tips	Bonus pack 10 x 50 mL cartridges, 30 blue mixing tips
<b>Futar®</b>	REF 1191211	REF 2827711
<b>Futar® Fast</b>	REF 1192611	REF 2827611
<b>Futar® Easy</b>	REF 1192311	-
<b>Futar® Easy Fast</b>	REF 1192411	-
<b>Futar® D</b>	REF 1193211	REF 2827811
<b>Futar® D Fast</b>	REF 1196111	REF 2827911
<b>Futar® D Slow</b>	REF 1195111	-



Made  
in  
Germany

Kettenbach GmbH & Co. KG  
Im Heerfeld 7  
35713 Eschenburg, Germany  
Phone +49 (0)2774 7050  
info@kettenbach.com  
www.kettenbach-dental.com



022744\_0523

  
**KETTENBACHDENTAL**  
Simply intelligent