




# COMFORTABLE UND PRECISE – **IDENTIVM®**

**VINYLSILOXANETHER® FOR PRECISION IMPRESSIONS**



**KETTENBACHDENTAL**  
Simply intelligent



# IDENTIUM® – THE UNIQUE MATERIAL COMBINATION!

## A special generation of precision impression material

Combines the advantages of two proven impression materials (A-silicone and polyether) – Vinylsiloxanether®.




### ONE-STEP IMPRESSION

- Double-mix impressions
- Monophase impressions




### INDICATIONS

- Crowns and bridges
- Inlays and onlays
- Functional impressions
- Implants
- Combined restorations (fixation impressions)
- Pick-up impressions



Double-mix impression




Implant impression



# IDENTIUM® – COMFORTABLE IN ALL DIMENSIONS

## Removal from the mouth - as light as a feather

Identium® vs. Impregum™ : 2.6 times lower withdrawal forces



Ease of removal from the patient's mouth avoids distortions and ensures a stress-free impression taking process.


**For your precision.**



**Advantages your patient will appreciate too**

## Fast final setting time

with sufficient processing time and short intraoral setting time



Sufficient intraoral working time makes it ideal for large cases.  
**For your precision.**

## Outstanding flowability



With the outstanding flowability of **Identium® Light**, you can reproduce even the narrowest gaps in great detail. **For your precision.**



- **Odorless and tasteless material** reduces gag reflex
- **Quick setting and easy removal from the mouth** for a relaxed and pleasant experience
- **Initially hydrophilic** – good performance even with residual moisture



# COMFORTABLE AND PRECISE IMPRESSIONS



EASY REMOVAL FROM THE MOUTH

SHORT INTRAORAL SETTING TIME

INITIALLY HYDROPHILIC


Available in viscosities  
**Heavy, Medium and Light**



Selectable setting times  
**regular** (4:30 min.)  
**fast** (3:30 min.)



Compatible with most common  
**mixing devices**



# IDENTIUM® – THE RANGE

HEAVY	Identium® Heavy		Identium® Heavy Fast		
	Refill pack	REF 14725 2 × 380 mL cartridges	REF 14727 2 × 380 mL cartridges		
MEDIUM		Identium® Medium		Identium® Medium Fast	
	Refill pack	REF 14717 2 × 380 mL cartridges	REF 14719 2 × 380 mL cartridges		
LIGHT		Identium® Light		Identium® Light Fast	
	Normal pack	REF 13701 2 × 50 mL cartridges, 8 gray mixing tips	REF 13711 2 × 50 mL cartridges, 8 gray mixing tips		
ACCESSORIES		Universal Adhesive	Dynamic Mixers	Mixing tips	Sympress Dispenser
	REF 14206  10 ml bottle	REF 17902  45 pieces white	REF 17226 / REF 17240  50 mixing tips gray / 100 mixing tips gray	REF 35911  Electrical dosing and mixing device	REF 17204  2 application syringes, 12 syringe tips



Kettenbach LP  
 16052 Beach Blvd. Ste. 221  
 Huntington Beach, CA 92647, USA  
 Phone: 877-532-2123  
 Fax: 888-892-9820  
 Website: [www.kettenbach-dental.us](http://www.kettenbach-dental.us)  
 Web shop: [www.kettenbachusa.com](http://www.kettenbachusa.com)  
 Custservice@kettenbachusa.com

  
**KETTENBACHDENTAL**  
 Simply intelligent