




TRUST IN PRECISION – **PANASIL®**

A-SILICONE FOR PRECISION IMPRESSION TAKING



KETTENBACHDENTAL

Simply intelligent

PRECISION WITHOUT COMPROMISES

Precision impression material (A-silicone)

In the **Panasil®** product family you find perfectly coordinated, innovative materials for all impression techniques and indications.

DOUBLE-MIX IMPRESSIONS



Panasil® tray Soft Heavy


PLEASANTLY HEAVY-BODIED



Panasil® initial contact Light

LOW VISCOSITY

INGENIOUSLY HYDROPHILIC



FOR EVERY TECHNIQUE AND EVERY INDICATION

TWO-STEP IMPRESSIONS

THE FAVOURITE



Panasil® tray Fast Heavy

SHORT SETTING TIME



Panasil® Putty Fast

**GOOD VALUE FOR
MONEY**



Panasil® binetics Putty Fast


**REAL PUTTY MATERIAL,
VERY EASY TO TRIM**



Panasil® initial contact X-Light

VERY LOW VISCOSITY

INGENIOUSLY HYDROPHILIC



PUTTY/WASH IMPRESSIONS



Panasil®
binetics Putty Soft


Panasil® Putty Soft

SMOOTH KNEADABLE, STABLE



Panasil® initial contact Light

**LOW VISCOSITY,
INGENIOUSLY HYDROPHILIC**



MONOPHASE IMPRESSIONS



Panasil® monophase Medium

HIGH FINAL HARDNESS



Panasil® monophase Medium

HIGH FINAL HARDNESS




OUTSTANDINGLY HYDROPHILIC FROM THE 1ST SECOND

Contact angle measurement


(parameter for the hydrophilicity of a material)

The flatter the drop, the more pronounced the hydrophilicity of the material.



Panasil® initial contact shows extremely fast spreading of the water drop and reaches the lowest contact angle after 3 seconds compared to the competitive materials.

EXCELLENT HYDROPHILICITY



- **Panasil® immediately overcomes any residual moisture** on direct contact with the tooth surface and ensures a precise result.



- **Panasil® ensures sharply defined preparation margins** even in **extreme situations**, e.g. in the case of difficult hemostasis.

PANASIL® – THE RANGE

PUTTY	Panasil® Putty		Panasil® Putty Soft	Panasil® Putty Fast
	Normal pack	REF 11101	REF 11121	REF 11141
		450 mL base paste, 450 mL catalyst paste, 2 dispensing scoops		
	Panasil® binetics Putty Soft		Panasil® binetics Putty Fast	
	Refill pack	REF 14703		REF 14701
		2 × 380 mL cartridges		

HEAVY	Panasil® tray Soft Heavy		Panasil® tray Fast Heavy	Panasil® tray Soft Heavy Fast
	Refill pack	REF 14707	REF 14705	–
		2 × 380 mL cartridges		
	Normal pack	REF 13541	REF 13551	REF 13561
		2 × 50 mL cartridges, 6 white mixing tips		

MEDIUM	Panasil® monophasic Medium			
	Refill pack	REF 14709		
		2 × 380 mL cartridges		
	Normal pack	REF 13501		
		2 × 50 mL cartridges, 6 white mixing tips		

LIGHT	Panasil® initial contact X-Light/ X-Light Fast		Panasil® initial contact Light/ Light Fast	Panasil® initial contact Regular
	Normal pack	REF 13401/ REF 13461	REF 13411 / REF 13471	REF 13431
		2 × 50 mL cartridges, 8 gray mixing tips		2 × 50 mL cartridges, 6 white mixing tips

Kettenbach LP
16052 Beach Blvd. Ste. 221
Huntington Beach, CA 92647, USA
Phone: 877-532-2123
Fax: 888-892-9820
Website: www.kettenbach-dental.us
Web shop: www.kettenbachusa.com
Custservice@kettenbachusa.com

022742_2525



KETTENBACHDENTAL
Simply intelligent