



MUCOPREN® SOFT

Prosthetic relining


Mucopren® Soft is a soft long-term relining material based on A silicone for **direct and indirect application**.



Special long-term adhesion:
Excellent adhesion to acrylic dentures through patented Mucopren® Adhesive



Unique sealing through its very smooth, hydrophobic silicone surface:
Inhibits bacterial colonization thus ensuring hygienic conditions on the prosthesis¹⁾



Simple treatment & fast processing

Faster setting times and shorter waiting times ensure a fast and smooth process in the practice and laboratory.

Supports the healing process of implants

Also suitable for use in implantology for trouble-free osseointegration.



High wearing comfort

The definitive resilient reline is completely odorless and tasteless and makes wearing the prosthesis as pleasant as possible.

¹⁾ Result from in-vitro studies in which it was shown that the growth of *Candida albicans* was prevented thanks to the smooth silicone surface



KETTENBACHDENTAL
Simply intelligent

HOLDS TO ITS PROMISES



FURTHER INDICATIONS



- For the creation of **permanently soft relinings**
- For relieving pressure on the surgical site after insertion of **implants**
- For balancing mucous membrane areas with different resilience
- For relieving sharp bony ridges
- For relining **interim prostheses**
- As **mucous membrane conditioner** in case of prosthetic stomatitis
- For creating and relining **obturators**

Benefits of using the new adhesive

» ONLY A SINGLE
APPLICATION
NECESSARY



CHAIRSIDE

3 MIN

The processing time is reduced to only **3 minutes** following mouth removal. No water bath necessary.



LABSIDE




3 MIN

The processing time is reduced to only **3 minutes** in the water bath or pressure pot.

MUCOPREN® SOFT – THE RANGE



Mucopren® Soft

Mucopren® Soft Basic Set	REF 28105	50 mL Mucopren® Soft, 50 mL silicone sealant, 10 mL adhesive, 6 blue mixing tips, 16 red mixing tips, 1 brush holder, 20 single-use brushes, 1 steel bur, accessories	
Mucopren® Soft Normal pack	REF 15687	2 x 50 mL Mucopren® Soft, 6 blue mixing tips	
Mucopren® Silicone sealant	REF 15686	50 mL Mucopren® Silicone sealant, 8 red mixing tips	
Mucopren® Adhesive	REF 14203	10 mL Mucopren® Adhesive	



Technical data

Product name	Mixing ratio	Working time at 73 °F ≤	Intraoral setting time at 95 °F ≥	Water bath	End of setting ≥	Shore hardness (approx.)
Mucopren® Soft – direct relining	1:1	2:15 min.	3:15 min.	Not applicable	5:30 min.*	A 28
Mucopren® Soft – indirect relining	1:1	2:15 min.	Not applicable	3:00 min.	> 3 min.	A 28
Mucopren® Soft silicone sealant	1:1	1:30 min.	Not applicable	Not applicable	6:30 min. at 23°C	A 33

*Total setting time (removal time from mouth) from start of mixing

NOTES

Kettenbach LP
 16052 Beach Blvd. Ste. 221
 Huntington Beach, CA 92647, USA
 Phone: 877-532-2123
 Fax: 888-892-9820
 Website: www.kettenbach-dental.us
 Web shop: www.kettenbachusa.com
 Custservice@kettenbachusa.com

Not for use, sale and distribution in Germany.



KETTENBACHDENTAL
Simply intelligent