



SOOO ADHERENT IS MUCOPREN® SOFT

Holds to the denture. And does what it promises.

Mucopren® Soft is a permanently soft durable relining material for removable dentures based on vinyl polysiloxane. The material can be applied directly or indirectly, and is also very useful for many other indications, like relieving the operation site after placing implants, for compensating areas of mucous membrane with varying degrees of resiliency and f.e. also for relieving sharp bony ridges.



User advantages

- Resistant to contamination*
- Can be precisely trimmed
- Excellent bond strength to the denture
- For direct and indirect application

Product properties

- Smooth, hydrophobic silicone surface*
- Silicone-based sealing
- High tear-resistance
- Patented adhesive

Benefits for the dentist

- Optimum requirements for hygienic denture conditions
- Long period in situ and comfort
- Modern material characteristics ensure easy, reliable handling
- Very efficient, value creation remains in the dental practice

For denture relining
without compromises!



Direct reline



Indirect reline



Can be easily trimmed



Easy to seal



* Concluded from in vitro studies which showed that the smooth silicone surface prevented *Candida albicans* growth.