

THE LUTING COMPOSITE FOR
ALL INDIRECT
RESTORATIONS



SOOO SMART
IS **VISALYS**[®] CEMCORE

CEMENTATION & CORE BUILD-UP



KETTENBACHDENTAL
Simply intelligent

VISALYS® CEMCORE – ENTIRELY CONVINCING

FOR ALL CASES

- **Suitable for permanent cementation** of crowns, bridges, inlays, onlays, veneers and adhesive bridges, root posts and indirect core build-ups
- **Suitable for all dental restoration materials** ((noble)metal, oxide ceramic, (zirconium oxide and aluminum oxide ceramics), silicate ceramic (feldspar and glass ceramic), hybrid ceramic, composite)

Conclusion by Prof. Dr. Claus-Peter Ernst

*[...] (One) impressive potential by **Visalys® CemCore** for a huge range of indications [...]*

Source: Ernst, Prof. Dr. Claus-Peter et al.: A universal adhesive cementation and core build-up material, in: Spitta: ZMK, Vol. 36, Issue 5 May 2020 [294–308]



HIGH COLOR STABILITY

Unlike other luting composites, **Visalys® CemCore** is hydrophobic. This minimizes the absorption of colors and discoloration of the cement joint.



	RelyX™ Universal	RelyX™ Unicem	RelyX™ Ultimate	Panavia™ F2.0	Panavia™ V5	Multilink® Automix	Variolink® Esthetic DC	Visalys® CemCore
Water								
Coffee								
Red wine								

Light-cured samples, stored 48 h in the corresponding solution at 37 °C.

SECURE ADHESION & LONG LIFE





The **high quality of the marginal gap** of **Visalys® CemCore** represents the high quality of the bond between the dental substance and the composite.

Source: Prof. Dr. Roland Frankenberger et al. (2019), *Marginal Quality of Ceramic Inlays luted with Visalys® CemCore*, University Marburg

No debonding and no fracture, even after **10 years** of simulated oral use.

Source: Prof. Dr. Dipl.-Ing. Martin Rosentritt (2020), *Fracture force after long-term storage and TCML*, University Regensburg



SMART HANDLING

Always the same procedure, regardless of the type of indirect restoration.



Easy to remove any excess

2)



**Now view the video
"Fixation of 16 veneers"**

1) Outstanding adhesion results with patented Active-Connect-Technology.
2) Copyright: Dr. med. dent. Frank-Michael Maier

VISALYS® CEMCORE – THE RANGE



Visalys® CemCore

Starter pack

4.5 g (2.5 mL) automix syringe,
2 mL Tooth Primer, 2 mL Restorative Primer,
blue mixing tips, blunt/tapered, 5 units each,
3 intraoral tips, 2 endo tips

Universal (A2/A3)

REF 13570

Normal pack

9 g (5 mL) automix syringe,
blue mixing tips, blunt/tapered, 10 units each,
6 intraoral tips, 4 endo tips

Universal (A2/A3)

REF 13572

Opaque

REF 13573

Translucent

REF 13574

Bleach

REF 13575

Dark (A4)

REF 13576

Visalys® Tooth Primer

4 ml bottle

REF 13580

Visalys® Restorative Primer

4 ml bottle

REF 13581

Visalys® CemCore Try In Paste

1.4 mL syringe,
5 application tips

Universal (A2/A3)

REF 13592

Opaque

REF 13593

Translucent

REF 13590

Bleach

REF 13591

Dark (A4)

REF 13594



NOTES

Kettenbach LP
16052 Beach Blvd. Ste. 221
Huntington Beach, CA 92647, USA
Phone: 877-532-2123
Fax: 888-892-9820
Website: www.kettenbach-dental.us
Web shop: www.kettenbachusa.com
Custservice@kettenbachusa.com

0210524_0724



KETTENBACHDENTAL
Simply intelligent