




SOOO STABLE
IS **VISALYS**[®] CORE

CORE BUILD-UP & POST CEMENTATION



KETTENBACHDENTAL

Simply intelligent

VISALYS® CORE – EASY TREATMENT

Dual-curing core build-up & root post cementation material

Visalys® Core – especially for difficult and freehand core build-ups.



EXCELLENT STABILITY

- with very good flow properties at the same time



UNIQUE ACTIVE-CONNECT-TECHNOLOGY

- for optimized adhesive bonding
- for use with common adhesives

DENTIN-LIKE GRINDABILITY




FREE OF BISPHENOL A



Distinguished by the experts “The Dental Advisor”


VISALYS® CORE – STABLE RESULTS

Adhesive bonding between Visalys® Core and dentin with different adhesives



Visalys® Core achieves the clinically relevant adhesive bond of approx. 10 MPa.

Average removal depth when grinding various core build-up materials



Visalys® Core can be precisely processed similar to dentin.

VISALYS® CORE – THE RANGE

Visalys® Core

Normal pack

White
Dentin

2 × 9 g (5 mL) automix syringe,
20 brown mixing tips
10 intraoral tips, 10 endo tips

Normal pack

White
Dentin

1 × 45 g (25 mL) cartridge,
20 yellow mixing tips,
20 intraoral tips



NOTES

Kettenbach LP
16052 Beach Blvd. Ste. 221
Huntington Beach, CA 92647, USA
Phone: 877-532-2123
Fax: 888-892-9820
Website: www.kettenbach-dental.us
Web shop: www.kettenbachusa.com
Custservice@kettenbachusa.com

0210518_0724



KETTENBACHDENTAL
Simply intelligent