Printing date 02/17/2023

Reviewed on 08/19/2021

## **1** Identification

- · Product identifier
- Trade name: Visalys Flow
- · Application of the substance / the mixture Dental filling material
- · Details of the supplier of the safety data sheet
- *Manufacturer/Supplier: Kettenbach GmbH & Co. KG Im Heerfeld 7 35713 Eschenburg / Germany*
- Information department: Regulatory Affairs Tel. +49 (0)2774/705-0
- Emergency telephone number: Regulatory Affairs Tel. +49 (0)2774/705-0

## 2 Hazard(s) identification

Health effects are based on the unmixed/uncured material. Under FDA regulations this product is classified as a medical device. Medical devices are exempt from the labelling provisions of GHS according to the Hazard Communication Standard (HCS).

### · Classification of the substance or mixture



Skin Sens. 1 H317 May cause an allergic skin reaction.

- · Label elements
- · Hazard pictograms



### · Hazard statements

- May cause an allergic skin reaction.
- · Precautionary statements
- *Wear protective gloves / protective clothing.*
- If on skin: Wash with plenty of water.
- If skin irritation or rash occurs: Get medical advice/attention.
- Dispose of contents/container in accordance with local/regional/national/international regulations.
- · Other hazards Use only for intended scope of application
- · Results of PBT and vPvB assessment
- · **PBT:** Not applicable.
- · vPvB: Not applicable.

## 3 Composition/information on ingredients

### · Chemical characterization: Mixtures

• Description: Mixture of the substances listed below with nonhazardous additions.

· Dangerous	components.
Dungerous	components

· Dangerous comp	onents:	
CAS: 6606-59-3	1,6-hexanediylbismethacrylate	5-10%
	𝔣 Skin Irrit. 2, H315; Eye Irrit. 2A, H319; STOT SE 3, H335	
CAS: 1565-94-2	BIS GMA	2.5-5%
	🚸 Skin Irrit. 2, H315; Eye Irrit. 2A, H319; Skin Sens. 1, H317	
CAS: 109-16-0	Triethylene glycol dimethacrylate	<i>≤</i> 2.5%
	🚸 Skin Sens. 1, H317	
CAS: 72869-86-4	Urethanedimethacrylate	<i>≤</i> 2.5%
	🚸 Skin Sens. 1B, H317	

(Contd. on page 2)

Printing date 02/17/2023

Reviewed on 08/19/2021

Trade name: Visalys Flow

(Contd. of page 1)

## 4 First-aid measures

- · Description of first aid measures
- General information: No special measures required.
- · After inhalation: Seek medical treatment in case of complaints.
- After skin contact: Immediately wash with water and soap and rinse thoroughly.
- After eye contact: Rinse opened eye for several minutes under running water. Then consult a doctor.
- · After swallowing: If symptoms persist consult doctor.
- Information for doctor:
- Most important symptoms and effects, both acute and delayed No further relevant information available.
- · Indication of any immediate medical attention and special treatment needed
- No further relevant information available.

## 5 Fire-fighting measures

- · Extinguishing media
- Suitable extinguishing agents: Use fire fighting measures that suit the environment.
- Special hazards arising from the substance or mixture No further relevant information available.
- · Advice for firefighters
- · Protective equipment: No special measures required.

#### 6 Accidental release measures

- · Personal precautions, protective equipment and emergency procedures Not required.
- · Environmental precautions: No special measures required.
- · Methods and material for containment and cleaning up: Pick up mechanically.
- · Reference to other sections

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

- See Section 13 for disposal information.
- · Protective Action Criteria for Chemicals

· PAC-1:		
CAS: 109-16-0	Triethylene glycol dimethacrylate	33 mg/m <sup>3</sup>
CAS: 72869-86-4	Urethanedimethacrylate	120 mg/m <sup>3</sup>
· PAC-2:		
CAS: 109-16-0	Triethylene glycol dimethacrylate	360 mg/m <sup>3</sup>
CAS: 72869-86-4	Urethanedimethacrylate	1,300 mg/m <sup>3</sup>
· PAC-3:		
CAS: 109-16-0	Triethylene glycol dimethacrylate	2,100 mg/m <sup>3</sup>
CAS: 72869-86-4	Urethanedimethacrylate	7,900 mg/m <sup>3</sup>

## 7 Handling and storage

· Handling:

### · Precautions for safe handling

For dental use only.

No special precautions are necessary if used correctly.

Observe the instructions for use.

• Information about protection against explosions and fires: No special measures required.

· Conditions for safe storage, including any incompatibilities

· Storage:

• Requirements to be met by storerooms and receptacles: No special requirements.

(Contd. on page 3)

<sup>–</sup> US

Printing date 02/17/2023

Reviewed on 08/19/2021

Trade name: Visalys Flow

- · Information about storage in one common storage facility: Not required.
- $\cdot$  Further information about storage conditions: See storage instructions on package
- $\cdot$  Specific end use(s) No further relevant information available.

## 8 Exposure controls/personal protection

• Additional information about design of technical systems: No further data; see item 7.

- · Control parameters
- · Components with limit values that require monitoring at the workplace:
- The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- Additional information: The lists that were valid during the creation were used as basis.
- · Exposure controls
- · Personal protective equipment:
- General protective and hygienic measures:
- The usual precautionary measures for handling chemicals should be followed.
- · Breathing equipment: Not required.
- · Protection of hands:



Protective gloves

#### · Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked prior to the application.

· Eye protection: Safety glasses

# 9 Physical and chemical properties

General Information		
Appearance:	X 7:	
Form:	Viscous	
Color:	According to product specification	
Odor:	Characteristic	
Odor threshold:	Not determined.	
pH-value:	Not determined.	
Change in condition		
Melting point/Melting range:	Undetermined.	
Boiling point/Boiling range:	Undetermined.	
Flash point:	Not applicable.	
Flammability (solid, gaseous):	Not applicable.	
Decomposition temperature:	Not determined.	
Auto igniting:	Product is not selfigniting.	
Danger of explosion:	Product does not present an explosion hazard.	
Density:	Not determined.	
Relative density	Not determined.	
Vapor density	Not determined.	
Evaporation rate	Not determined.	

Printing date 02/17/2023

Reviewed on 08/19/2021

Trade name: Visalys Flow

		(Contd. of page
· Solubility in / Miscibility with		
Water:	Not miscible or difficult to mix.	
· Partition coefficient (n-octan	ol/water): Not determined.	
· Viscosity:		
Dynamic:	Not determined.	
Kinematic:	Not determined.	
• Other information	No further relevant information available.	

## 10 Stability and reactivity

· Reactivity No further relevant information available.

· Chemical stability

• Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.

- · Possibility of hazardous reactions No dangerous reactions known.
- · Conditions to avoid No further relevant information available.
- · *Incompatible materials:* No further relevant information available.
- Hazardous decomposition products: No dangerous decomposition products known.

### 11 Toxicological information

· Information on toxicological effects

- · Acute toxicity:
- · Primary irritant effect:
- on the skin:

*Testing according to ISO 10993-1/10993-10 showed no skin irritation for this product when fully cured* • *on the eye:* No irritating effect.

· Sensitization: Sensitization possible through skin contact.

· Additional toxicological information:

The product is not subject to classification according to internally approved calculation methods for preparations:

When used and handled according to specifications, the product does not have any harmful effects according to our experience and the information provided to us.

#### · Carcinogenic categories

· IARC (International Agency for Research on Cancer)		
CAS: 7631-86-9	silicon dioxide, chemically prepared	3
CAS: 128-37-0	Butylated hydroxytoluene	3
· NTP (National T	Foxicology Program)	
None of the ingre	edients is listed.	
· OSHA-Ca (Occupational Safety & Health Administration)		

None of the ingredients is listed.

## 12 Ecological information

· Toxicity

- Aquatic toxicity: No further relevant information available.
- · Persistence and degradability No further relevant information available.
- $\cdot$  Behavior in environmental systems:
- · Bioaccumulative potential No further relevant information available.
- *Mobility in soil* No further relevant information available.

(Contd. on page 5)

US

Printing date 02/17/2023

Reviewed on 08/19/2021

#### Trade name: Visalys Flow

- · Additional ecological information:
- · General notes:
- Water hazard class 1 (Self-assessment): slightly hazardous for water Do not allow product to reach ground water, water course or sewage system.
- · Results of PBT and vPvB assessment
- · PBT: Not applicable.
- · vPvB: Not applicable.
- · Other adverse effects No further relevant information available.

## 13 Disposal considerations

- · Waste treatment methods
- · Recommendation:

Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

· Uncleaned packagings:

· Recommendation: Disposal must be made according to official regulations.

4 Transport information	
· UN-Number · DOT, ADR, ADN, IMDG, IATA	Void
· UN proper shipping name · DOT, ADR, ADN, IMDG, IATA	- Void
· Transport hazard class(es)	-
· DOT, ADR, ADN, IMDG, IATA · Class	Void
· Packing group · DOT, ADR, IMDG, IATA	- Void
• Environmental hazards: • Marine pollutant:	No
· Special precautions for user	Not applicable.
• Transport in bulk according to Annex II MARPOL73/78 and the IBC Code	of Not applicable.
· UN "Model Regulation":	Void

## 15 Regulatory information

- · Safety, health and environmental regulations/legislation specific for the substance or mixture · Sara
- · Section 355 (extremely hazardous substances):

None of the ingredients are listed.

### · Section 313 (Specific toxic chemical listings):

None of the ingredients are listed.

#### • TSCA (Toxic Substances Control Act):

One or more substances contained in this material are not listed on the TSCA Inventory. Contact Kettenbach for more information.

### · Hazardous Air Pollutants

None of the ingredients is listed.

(Contd. on page 6)

(Contd. of page 4)

<sup>•</sup> US

(Contd. of page 5)

## Safety Data Sheet acc. to OSHA HCS

Printing date 02/17/2023

Reviewed on 08/19/2021

Trade name: Visalys Flow

· Proposition 65

 $\cdot$  Chemicals known to cause cancer:

None of the ingredients are listed.

· Chemicals known to cause reproductive toxicity for females:

None of the ingredients are listed

· Chemicals known to cause reproductive toxicity for males:

None of the ingredients are listed

· Chemicals known to cause developmental toxicity:

None of the ingredients are listed

· Cancerogenity categories

· EPA (Environmental Protection Agency)

None of the ingredients are listed

· TLV (Threshold Limit Value established by ACGIH)

Contact Kettenbach for more information.

· NIOSH-Ca (National Institute for Occupational Safety and Health)

None of the ingredients are listed

### 16 Other information

The information contained herein is based on our present knowledge and is believed to be correct at the date issued. The information is subject to change without notice. However, all information, statements, recommendations or suggestions shall not constitute any guarantee, representation or warranty, expressed or implied.

Final determination whether the product is fit and suitable for a particular purpose and the user's method of use or application is the sole responsibility of the user.

· Department issuing SDS: Regulatory Affairs

- · Contact: Regulatory Affairs Manager Tel. +49 (0)2774/705-0
- · Date of preparation / last revision 08/19/2021 / 5
- Abbreviations and acronyms:

Skin Irrit. 2: Skin corrosion/irritation – Category 2

Eye Irrit. 2A: Serious eye damage/eye irritation – Category 2A Skin Sens. 1: Skin sensitisation – Category 1

Skin Sens. 1B: Skin sensitisation – Category 1B

STOT SE 3: Specific target organ toxicity (single exposure) - Category 3

us -